




# Olympus FVMPE-RS Hybrid Multiphoton System Standard Operation Protocol Basic Operation




## Turn on system



Please sign on the log sheet before switching on system.


- Turn on Chameleon laser by turning the Laser Key 90° clockwise to horizontal position (if only one laser needs to be used, please only turn on laser 1)
- Turn on computer and log in with “User” ①
- Turn on power supply Unit FV30-PSU ② (please wait for about 10 sec before next step)
- Turn on power switch of Main scanner ③ (please wait for about 10 sec before next step)
- Turn on power switch of SIM scanner ④
- Turn on fluorescent lamp illuminator ⑤ (the switch is at the left back of the box, shown by →), wait for 3sec, press the “ON/OFF button at the front panel
- Turn on microscope power ⑥
- Turn on the touch panel controller ⑦ (the switch is at the right back of the box, shown by →)

## Setting the microscope


- Double click the “FV30S-SW” icon  on the desktop to initialize the software, it may take 1-2min for software and hardware connection.
- Keep pressing  on touch screen to move the objective up
- Select the objective lens manually and engage it in the light path 
- Put your sample on the stage
- Select **Ocular** to set the microscope for direct observation via the eyepieces
- Click on **EPI** for fluorescence observation
- Select the filter for fluorescence observation in the drop down list of **Cube Turret #2**.




- Put water on the sample if water objective will be used.
- Keep pressing  on touch screen to move the objective close to the sample
- Click **Open** in **EPI shutter** to open shutter if necessary 
- Focus the sample with Coarse/fine focus knob
- Locate the view of interest through stage control knob




## Configuring the beam path and laser


- In **IR Laser Setting** Tool Window, set the wavelength of the IR laser to be used.
- Press the button  of **Fine adjustment**, and press the button **Start active alignment** which will align the laser to optimal position.




- Click on the **Lightpath** button to open the **Lightpath** window.



- Select coupling mirror in **Coupling mirror of IR laser** tool window depends on the wavelength to be used.
  - Chose "**Glass**" if only laser 1 is in use.
  - Chose "**LCDM 690-870**" for two laser lines with Laser 1  $\geq 920\text{nm}$  and Laser 2  $\leq 870\text{nm}$
  - Chose "**LCDM 1030-1300**" for two laser lines with Laser 1  $\leq 940\text{nm}$  and Laser 2  $\geq 1030$



- Check NDD filters in **NDD** panel of Lightpath





|        |   |
|--------|---|
| SDM    |   |
| SDM570 | >570nm pass to Ch1/2, <570nm reflect to Ch3/4 |
| Empty  | All pass to Ch1/2                             |
| Mirror | All reflect to Ch3/4                          |

**NDD filter cube list at FCF:**



| FCUBE           | Dichroic mirror | BF for Ch1/Ch3      | BF for Ch2/Ch4   |
|-----------------|-----------------|---------------------|------------------|
| Blue/Green      | DM485           | G 495-540 (Green)   | B 410-460 (Blue) |
| Red/FarRed      | DM650           | FR 660-750 (FarRed) | R 575-645 (Red)  |
| GFP/mCherry     | DM560           | mCherry LP590       | GFP 525/50       |
| mCherry/Qdot655 | DM630           | Qdot 655/40         | mCherry 600/50   |
| SHG/YFP         | DM495           | YFP 535/50          | SHG 460/50       |

*Please be noted that FCUBE filter might be changed based on user's application. Please contact FCF staff if the FCUBE filter in use does not fit your application.*

- Press the **Dye & Detector Select** button in **PMT Setting** tool window. The **Dye & Detector Select** dialog box appears
- Select the detector to be assigned to the observation channel in Detector list
- Press the **Add** button. The observation channel number and the abbreviation (physical channel name) of the assigned detector are displayed in the observation channel list.
- After setting all channels, press the button **OK**.




- If sequential scan needed for multi-channel acquisition, click **Line** or **Frame** in **PMT Setting** tool window





- **Simultaneous scan**
  - Use for one or two excitation wavelength
  - Advantage: faster image acquisition
  - Disadvantage: cross talk might be happened between channels

- **Sequential scan**
  - Only applicable for two different laser wavelength
  - Advantage: when one group is active, only one detector and one laser line is switched on. This may reduce cross talk
  - Disadvantage: slower image acquisition

**Reuse previous acquisition settings**

- Select **Open - Image** in the **File** menu to display the dialog box, and select the image to view.
- Click **Show/Hide** and select **Property** to load information of acquisition settings
- Click  to reload previous setting




**Adjusting the live image**

- Click on the **LSM** button to switch to scanning mode, then click **LSM imaging** for scan settings



- In **PMT** window, select the laser, set an initial value for laser intensity and HV value  
*It is suggested to set the laser intensity <5% and HV <800*
- Press the **Live** button in the **Live** window. (The button view is changed to **Stop**.) The live image is displayed in display area.
- In **PMT** Window, fine tune Laser Intensity (%), Sensitivity (HV), Gain and Offset to optimize the image



- Click on **Hi-Lo** button to display saturation pixels. Red = Saturation (maximum). Blue = Zero (minimum).
- Adjust **laser intensity (%)** and **HV** to obtain few red pixels shown on image. Increase **Offset** to introduce blue pixels in background area.
- Click **Stop** to complete the preview of the sample.



\* User is advised **not to use too extreme settings on laser intensity and HV** (as **too high laser power will induce bleaching** onto the sample; **sample burning** might be happened if the laser power is too high; and too high HV will cause significant noise and shorten the lifetime of the PMT). In any given cases, users **should keep the HV < 800**.



## Acquiring an image

- In **LSM Imaging** Tool Window, select the scanner **Type** and **Mode**.



|             |   |
|-------------|---|
| <b>Type</b> | "Galvano" (High definition scan) or "Resonant" (High speed scan)          |
| <b>Mode</b> | "OneWay" (Scan in one direction) or "Roundtrip" (Scan in both directions) |

**Galvano** and **OneWay** are suggested in regular acquisition.


- Use the slider in the **Speed** to adjust the scan speed. Slower speed increases signal/noise ratio of the image.
- Select the **Scan Size** as predefined number of pixels.
- Select the number of **average** in **PMT** tool window. Averaging improves the image by increasing the signal-to-noise ratio.
- In "Acquire" tool window, press the  button and select your folder in **F** drive to save the images. The acquired images are saved automatically.
- Press  to start acquiring the image.









## Scanning Z Series

- Select **ON** for **Z** in **Series** Tool Window.
- Select "**Start/End**" in **Motor** on **Z Section**.
- Click  button to have continuously scanning.
- Press the  button in **Origin** at the Z position which is the reference position to acquire the images.  
*\*The distance from the Z position registered in [Origin] to the current Z position is displayed in [Current].  
If you register the top surface of the specimen in [Origin], the depth from the top surface to the observation position is always displayed in [Current], which is very useful.*




- Turn the focusing knob gently to locate the start position of the specimen to **Register** start. Then locate the end position of specimen to **Register** end.  
*Press the **START<=>END** button to reverse the start position and the end position if necessary.*


- Click  to complete the preview of the sample.
- Set "**Step Size**" for Z-stack scanning. *This digit suggests the distance in  $\mu\text{m}$  between each slide.*

- Click on the  button to start the recording of the Z Series.
- After the image is acquired, pressing the  button allows you to perform the image acquisition repetitively under the same condition. Press the  button to complete the image acquisition.




## Setting Time Series

- Select **ON** for **Time** in **Series** Tool Window.
- Set the interval by input in how many seconds in **Interval**.  
*\*If you attempt to set the time shorter than the time displayed in [Scan] (time taken for acquiring 1 cycle image) in [Interval], "FreeRun" appears. In this case, the interval to acquire the image is the time displayed in [Scan].*
- Set the number of image acquisitions in **Cycle**.  
*The time from the start of the T series image acquisition to the end of the T series image acquisition is calculated and displayed in [LSM Total]. Set [Cycle] so that [LSM Total] becomes an appropriate time*
- Click on the  button to start the recording of the Z Series.
- After the image is acquired, pressing the  button allows you to perform the image acquisition repetitively under the same condition. Press  the button to complete the image acquisition.



## Turn off system

**Please check if the equipment will be used by other users. Please switch off system if no one books equipment over two sessions (1h) after you.**

1. Soak the water objective in MilliQ water for 5min
2. Turn off fluorescent lamp illuminator by pressing "ON/OFF" button for **3 seconds**, the illuminator will show "300" to count down for **5min** for cooling down the lamp.
3. Turn off Chameleon laser by turning the laser key 90° anticlockwise to vertical position
4. Turn off power switch of SIM scanner ④
5. Turn off power switch of Main scanner ③
6. Turn off power supply Unit FV30-PSU ②
7. Exit "FV30S-SW" software and upload your data to server
8. Turn off computer ①
9. When the lamp illuminator **5min** count-down is finish, turn off the power switch of lamp illuminator ⑤ (the switch is at the left back of the box, shown by →)
10. Keep pressing  on touch screen to move the objective to upper limit, then clean up with lens paper
11. Turn off the touch panel controller ⑦ (the switch is at the right back of the box, shown by →)
12. Turn off microscope power ⑥

**\*Users are reminded to sign on the log sheet before departure.**